Instructions For Use

[Warning]

Before use, please contact your ophthalmologist. If the product has defect or broken during use, contact immediately your ophthalmologist for examination. Broken lens can bring ocular hazard e.g corneal erosion. However, proper use can bring following hazards inevitably.

- Possible decrease of corneal endothelial cell.
- Using contact lens interferes the oxygen supply, this can increase the possibility of ocular hazard e.g corneal neovascularization, corneal edema, keratitis and giant papillary hyperplasia.
- Product life time can be shortened depend on using condition. Change your contact lens periodically or discard when its lifetime is over.
- There are individual differences on usage time. If you feel discomfort, remove your contact lens. Extended use can decrease the supply of oxygen to corneal and bring the possibility of ocular hazards. Please use in proper time.
- Read carefully the attached instructions and follow the instructions for use, preservation and maintenance.

If you feel any problems, contact immediately you ophthalmologist for examination.

DO NOT USE your contact lens beyond the instructed purpose.

[Contraindications]

[Medical Contraindications]

- Acute and subacute inflammation of anterior
- Eye infection
- Uveitis
- Decreased corneal sensation
- Dry eye and lacrimal disorders on lens wearing
- Abnormal eyelid
- Allergy affecting lens wearing.

[General contraindications]

- Patients who cannot follow the medical instructions.
- Patients who cannot use in proper ways
- Patients who cannot get periodical inspection.
- Patients who cannot perform health control for lens use.
- Patients who are extremely nervous.

[Environmental contraindications]

- Patients who are in dry condition always.
- Patients whose eyes are exposed to dust and chemicals

[How to use]

This is daily wear contact lens. You must remove the lens before sleep. Don't forget the lens care when remove.

[Before use]

1. Check your eyes for gums and congestions with mirror.
2. Keep your nails short and round not to damage your eyes and contact lens.
3. Make sure that wash and rinse your hands.

[Checking for lens, package and container]

Before open, check the status of lens, package and container whether they have any defects. If you see any defects, don't open the container and contact the seller.

[Open]

1. Check the one lens is locating in the middle of the container, tear and status of color.
2. Turn over the cover, you can find arrow on the aluminum lid. Remove the lid twisting it in one direction. Be careful not to hurt your finger.
3. Remove the rubber cover and slightly shake the container to float the lens.
4. Poor it to lens case with preservative solution entirely.
5. If you failed to remove the lid, don't try excessively. Use the tools to open it.
6. Aged or infant who is not available to use the lens, do not try it excessively.

Poor it to lens case with preservative solution entirely.

4. If you failed to remove the lid, don't try excessively. Use the tools to open it.

[Confirmation]

1. Check the deformation, discoloration, debris, contamination, scratch and damages. If you are using lens with those defects, it can bring the possibility of ocular hazard. So DO NOT USE.
2. Put your lens on your finger and check in and outside. If you use the lens with inside out, you can feel discomfort and it easily come off.

[How to wear]

1. Put your lens on index finger of your dominant hand and check the in and outside of the lens. Moisture removal helps easy wearing.
2. Grab and fall down your lower eyelid with middle finger of same hand.
3. Lift up the upper eyelid with other hand's middle finger to raise the eye rims.
4. Put the lens on pupil with gazing your both eyes to the front.
5. Check the right location of the lens on pupil and slowly move back your finger.
6. Slowly blink your eyes for several times and check the lens was fitted rightly.

[How to remove]

1. Clean your hands before handling. Check the location of the lens with mirror and gaze the upper direction to make pupil go up.
2. Fall down the lower eyelid with middle finger of your dominant hand.
3. Grab and remove the lens with thumb and index finger of dominant hand. Be careful not to damage on to your eyeball and lens with nails.

[How to care]

- Daily chemical disinfects required. Do not use hot water.
- Please follow the attached instructions and marked information on the cover. Use this product in proper way.
- For treatment and preservation, use the solutions for soft contact lens. DO NOT tap water.

[Warnings on usage]

[Wearing warnings]

1. On wearing, follow the medical instructions.
2. If you feel pain, haziness and blurred focus during use, remove the lens immediately and contact your ophthalmologists.
3. Follow the medical instructions when need to use eye drops during wear.
4. To prevent from potential damage to corneal and lens, do not rub your eyes during wear.
5. During wear, watching TV and monitor for long time can bring the dry eye so blink your eyes repeatedly.
6. People who have allergic disease must contact ophthalmologist and follow instructions for proper and safe use before wear.
7. This contact lens has metallic ingredient pigment. If you are allergic to metal, contact your ophthalmologist. Remove the lens before MRI scanning.
8. Colored section partially blocks the light and it can make slight invisibility. This will be alleviating as continual use. If this keep being nervous, stop using and contact your ophthalmologist.
9. Until you get used to the lens, don't wear the lens during driving or operations of machinery. If you feel any problems during the task, make sure the safety and stop the operation.
10. Changes on the body like disease, pregnancy, childbirth and internal use of medicine, environment can affect on usage of lens. Please follow the medical instructions.
11. If you have any senses of congruity on your eye or feel bad, do not use the lens.
12. Aged or infant who is not available to use the lens, do not try it excessively.

Contact the ophthalmologist and follow the medical instruction.

13. Do not use other's lens. Do not rent your lens to anybody.

14. Do not lick the lens. Do not put the lens into your mouth.

15. Close your eyes when you use hairspray.

[Warnings on management]

1. Do not fold, bend the lens. And do not grab it with fingernail. If it's damaged, do not use.

2. Do not freeze the lens and its solution. Once the lens is frozen, its form is corrupted even if it melts down again. So do not use.

3. Do not use if debris, scars, transformation and discoloration are seen. 4. Do not let your lens dried. Once dried, it is not usable. Preserve it with solution.

Usage of the lens brings subjective symptoms sometimes as below. This will be alleviating as continual use. But if you feel bigger pain or these keep for long time, stop using and contact the ophthalmologist.

Symptom	Possible cause	How to cope
Sensation & Pain	Inside-out wear	Remove, clean, check and re wear the lens. Still have a problem, change it with new one.
	Debris	
	Contamination	
	Scars, damages	
	Ocular problem	Get the medical inspection and follow the instruction.
Ocular congestion	Long time use	Stop using, wait and see the eye condition. Rein on use of lens.
	Poor sleep, feeling bad	
	Contamination	Wash and disinfect the lens. If there is no progress, change it with new one.
	Dry eye	Get the medical inspection and follow the instruction.
	Allergy	
Ocular problem		
Itch on eyeballs	Contamination	Wash and disinfect the lens. If there is no progress, change it with new one.
	Allergy	Get the medical inspection and follow the instruction.
	Ocular problem	
Dry eye	Dry lens	Follow the medical instruction and use eye drops. Blink your eyes repeatedly.
	Dry air	
	Watching displaying monitors for long time	
	Dry eye	Get the medical inspection and follow the instruction.
	Medicine use	
Blurred View	Contamination	Wash and disinfect the lens. If there is no progress, change it with new one.
	Dry lens	Use eye drops.
	Feeling bad	Stop using, wait and see the eye condition. Rein on use of lens.
	Long time use	
	Ocular problem	Get the medical inspection and follow the instruction.

【 How to store 】

Keep away from a ray of sunlight. Store it at normal temperature.

【Maintenance】

Chemical disinfection is required.

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